



Switches

# Photonic Crystal Switch Rack




## Measurement Range

390 – 2000 nm

## Technical Specs

Available channels	4, 8
Guaranteed transmission	15 % @ 633 nm
Return loss	< 40 dB
Crosstalk	< -55 dB
Switching time	typ. 12 ms
Supply voltage	powered by wavelength meter (WR Series)
Operating temperature	0 – 60° C
Storage temperature	-40° – 80° C
Operating Relative Humidity	35 – 55 %
Dimensions	360 mm × 483 mm × 43.5 mm
Weight	4.0 kg
Standard interface	TTL, USB for service only





Switches

# Photonic Crystal Switch Rack



HighFinesse  
The Standard of Accuracy

## Further Information

For further technical information, application examples, diagrams and for customisation of photonic crystal switch racks please contact:

HighFinesse Service

[service@highfinesse.de](mailto:service@highfinesse.de)



HighFinesse GmbH  
Neckarsulmer Straße 5  
72072 Tübingen, Germany



T +49 (0) 7071 - 53 918 0  
F +49 (0) 7071 - 53 918 99  
M [info@highfinesse.com](mailto:info@highfinesse.com)



Additional information  
and distributors:  
[www.highfinesse.com](http://www.highfinesse.com)